




Medline International France SAS

5 rue Charles Lindbergh  
44110 Châteaubriant

[www.medline.com/en/](http://www.medline.com/en/)



OPS™ Ultimate

Drape sheet, three quarter

**Reference(s)**

**TB9349CEB**

**Description**

A sterile wrap 56x56 cm contains:  
- 1 tri-laminate three quarter drape, 145x193 cm

**Main materials**

Absorbent and impervious tri-laminated material (polyethylene + polypropylene)

**Pack factor**

16 units per carton

**Packaging**

Vented sterilization pouch  
1 polyliner contains all the units in a carton box

**Manufacturer**

Medline International France SAS  
EN ISO 13485

**Regulation**

In compliance with Medical Devices Directive 93/42/EEC  
MDD Class: I sterile  
HP: High performance according to EN 13795 standard



Latex free



0086



STERILE|EO



© 2014 Medline Industries Inc. or one of its subsidiaries.  
All rights reserved.

LIT-TB9349CEB-EN03

01.01.INT.12

FINT.441 Revision n°02 Effective date : 08 Jan 14