




Medline International France SAS
5 rue Charles Lindbergh
44110 Châteaubriant
www.medline.com/en/



OPS™ Ultimate

Drape sheet, three quarter

Reference(s)	TB9349CEB
Description	A sterile wrap 56x56 cm contains: - 1 tri-laminate three quarter drape, 145x193 cm
Main materials	Absorbent and impervious tri-laminated material (polyethylene + polypropylene)
Pack factor	16 units per carton
Packaging	Vented sterilization pouch 1 polyliner contains all the units in a carton box
Manufacturer	Medline International France SAS EN ISO 13485
Regulation	In compliance with Medical Devices Directive 93/42/EEC MDD Class: I sterile HP: High performance according to EN 13795 standard



© 2014 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-TB9349CEB-EN03

01.01.INT.12