




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



OPS™ Ultimate

Drape sheet, three quarter

Reference(s)

TB9349CEB

Description

A sterile wrap 56x56 cm contains:
- 1 tri-laminate three quarter drape, 145x193 cm

Main materials

Absorbent and impervious tri-laminated material (polyethylene + polypropylene)

Pack factor

16 units per carton

Packaging

Vented sterilization pouch
1 polyliner contains all the units in a carton box

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile
HP: High performance according to EN 13795 standard



Latex free



0086



STERILE EO



© 2014 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-TB9349CEB-EN03

01.01.INT.12

FINT.441 Revision n°02 Effective date : 08 Jan 14