




Medline International France SAS

5 rue Charles Lindbergh  
44110 Châteaubriant

[www.medline.com/en/](http://www.medline.com/en/)



## OPS™ Essential

## Head Drape

### Reference(s)

**ES29483CE**

### Description

A SMS wrap 66x66 cm contains:

- 1 bi-laminate head drape 240x150cm, 67x25cm reinforcement, 98 cm long adhesive, 1 cord holding tab

### Main materials

Absorbent and impervious bi-laminated material (polyethylene + polypropylene)  
Absorbent Spunbond polypropylene reinforcement

### Pack factor

38 units per carton

### Packaging

Paper and film peel pouch  
1 polyliner contains all the units in a carton box

### Manufacturer

Medline International France SAS  
EN ISO 13485

### Regulation

In compliance with Medical Devices Directive 93/42/EEC  
MDD Class: I sterile  
HP: High performance according to EN 13795 standard



Latex free



0086



© 2019 Medline Industries Inc. or one of its subsidiaries.  
All rights reserved.

TDS\_GenSurgDrape\_ES29483CE\_EN02

01.INT.03.30