




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



OPS™ Essential

Head Drape

Reference(s)

ES29483CE

Description

A SMS wrap 66x66 cm contains:

- 1 bi-laminate head drape 240x150cm, 67x25cm reinforcement, 98 cm long adhesive, 1 cord holding tab

Main materials

Absorbent and impervious bi-laminated material (polyethylene + polypropylene)
Absorbent Spunbond polypropylene reinforcement

Pack factor

38 units per carton

Packaging

Paper and film peel pouch
1 polyliner contains all the units in a carton box

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile
HP: High performance according to EN 13795 standard



Latex free



0086



© 2019 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS_GenSurgDrape_ES29483CE_EN02

01.INT.03.30