




Medline Industries Ltd
3rd Floor, Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL, United Kingdom

Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk



Technical Sheet

Reference(s)

Description

Main materials

Pack factor

Sterilization

Packaging

Legal Manufacturer

Regulation

ORIS Insufflation tubing Set

ORISISO, ORIS15, ORISLL, ORISPUSH

Medline ORIS insufflation tubing set is used to filter and insufflate clean CO₂ in order to establish pneumoperitoneum prior to laparoscopy procedures.

Insufflation tubing set is a single use device that comprises a flexible tubing and a filtering device.

The filtering device of Medline ORIS insufflation tubing set with >99,9% of Bacterial filtration efficiency (BFE) prevents from leaks, which minimize the risks of contamination of the patients and infections.

Item Code	Description	Length	Inner diameter
ORISISO	OR Insufflation Set with ISO 15M/22F connector	3 meters	7 millimeters
ORIS15	OR Insufflation Set with 15MM connector		
ORISLL	OR Insufflation Set with Luer-Lock connector		
ORISPUSH	OR Insufflation Set with female connector		

PVC, DEHP free, Latex free

25 units per case

Sterilized by Ethylene oxide

Double package

Advanced Medical Design Co., Ltd.,
4F~5F, No. 29, Wuquan 5th Rd., Wugu Dist.,
24888 New Taipei City, Taiwan.

In compliance with Medical Devices Directive 93/42/EEC
MD Class: IIa sterile

- Please refer to complete labeling and instruction for use (if applicable) before using the product.



© 2022 Medline Industries LP or one of its subsidiaries. All rights reserved.
TDS_ORIS_EN02