




Medline International France SAS
5 rue Charles Lindbergh
44110 Châteaubriant
www.medline.com/en/



OPS™ Ultimate

Drape sheet, half

Reference(s)	TB9358CE
Description	A sterile wrap 56x56 cm contains: - 1 tri-laminate drape, 112x145 cm
Main materials	Absorbent and impervious tri-laminated material (polyethylene + polypropylene)
Pack factor	20 units per carton
Packaging	Paper and film peel pouch 1 polyliner contains all the units in a carton box
Manufacturer	Medline International France SAS EN ISO 13485
Regulation	In compliance with Medical Devices Directive 93/42/EEC MDD Class: I sterile HP: High performance according to EN 13795 standard



© 2014 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-9358CE-EN02

01.01.INT.12