
Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



OPS™ Essential

Head Drape

Reference(s)	**ES29483CE**
Description	A SMS wrap 66x66 cm contains: - 1 bi-laminate head drape 240x150cm, 67x25cm reinforcement, 98 cm long adhesive, 1 cord holding tab
Main materials	Absorbent and impervious bi-laminated material (polyethylene + polypropylene) Absorbent Spunbond polypropylene reinforcement
Pack factor	38 units per carton
Packaging	Paper and film peel pouch 1 polyliner contains all the units in a carton box
Manufacturer	Medline International France SAS EN ISO 13485
Regulation	In compliance with Medical Devices Directive 93/42/EEC MDD Class: I sterile HP: High performance according to EN 13795 standard

Latex free

0086

STERILE EO

HP EN13795

© 2019 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.